

## *SO-LOW FLAMMABLE MATERIAL STORAGE FREEZER MODEL DHH20-14SDFMS MANUAL DEFROST*

*Meets NFPA and OSHA guidelines 45, and 70 for flammable material storage in general laboratory work areas*

### **14 Cubic Ft. Upright -25°C / Solid Door**

#### **Specifications:**

* Capacity	14 cubic ft. / 396 liters
* Temperature	-15°C to -25°C
* Dimensions	
* Exterior w x d x h	23.75" x 24.75" x 73"
*Interior w x d x h	19" x 17.5" x 63.75"

#### **Features:**

*Thermostat	*Microprocessor temperature controller
	*Remote alarm contacts
	*Power failure alarm
	*Audible and visual high/low temperature alarms
	*Alarm validation (audible, visual, and remote alarm)
	*Sample simulated digital temperature display
*Access Port	3/8" Access port
*Compressor	1/4 Hp Compressor
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	6 Fixed
*Insulation:	High density urethane foam
*Refrigerant:	R-600A
*Air Circulation	Cold Wall
*Defrost	Manual Defrost
*Voltage	115-60-1
*Running Amps	3.5
*Supply Plug	5-15 NEMA
*Shipping weight	195 lbs.
*Approvals	ETL Listed
*Warranty-(Domestic)	Two year parts and labor.

