

MODEL DHH20-14SDFX

EXPLOSION PROOF FREEZER -15°C TO -25°C

Capacity [Cu. Ft.] 14 Cu. Ft.
Capacity [Liters] 396 Liters
Temperature Range [°C] -15°C to -25°C
Temperature Range [°F] 5°F to -13°F

Dimensions Width x Depth x Height
 Exterior [Inches] 23.5" x 24.5" x 73"
 Interior [Inches] 19" x 17.5" x 63.75"
 Exterior [Centimeters] 60 x 62 x 185
 Interior [Centimeters] 48 x 44 x 162

Compressor Hermetic Compressor
Air Circulation Gravity Flow
Defrost Method Manual
Refrigerant R290
Exterior Painted Steel
Interior Painted Steel
Insulation High Density Urethane
Lockable Door 1 Door, Keyed Lock
Door Material Insulated Steel
Leveling Legs 4
Storage 6 Fixed Shelves
Shipping Weight 195 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 3.5A
Electrical Plug Not Supplied (Hardwire)
Certification(s) ETL Listed
 Energy Star Rated

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp
Min / Max Temperature Display and Reset
Alarm Relay Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» **Meets NFPA / OSHA Guidelines for Hazardous Locations Class 1, Division II, Group C & D**

Warranty 2 Years Parts & Labor.