

MODEL DHH20-17MDP

MANUAL DEFROST FREEZER -15°C TO -25°C

Capacity [Cu. Ft.] 17 Cu. Ft.
Capacity [Liters] 481 Liters
Temperature Range [°C] -15°C to -25°C
Temperature Range [°F] 5°F to -13°F

Dimensions Width x Depth x Height
 Exterior [Inches] 30.75" x 28.63" x 67.75"
 Interior [Inches] 26.13" x 22.5" x 58.25"
 Exterior [Centimeters] 78 x 72 x 172
 Interior [Centimeters] 66 x 57 x 147

Compressor Hermetic Compressor
 Air Circulation Gravity Flow
 Defrost Method Manual
 Refrigerant R290
 Exterior Powder Coating
 Interior Molded ABS
 Insulation Polyurethane Foam
 Lockable Door Keyed
 Porthole 3/4" Diameter
 Leveling Legs 4
 Storage 6 Fixed Shelves
 Shipping Weight 235 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
 Amperage Line | Running 15 Amp Dedicated | 5.5 A
 NEMA Configuration NEMA 5-15
 Certification(s) ETL Listed

Controller

Display Type Microprocessor, Digital
 Alarms High/Low Temp
 Alarm Relay Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 2 Year Parts & Labor