

MODEL DHH20-17MDP

MANUAL DEFROST FREEZER -15°C TO -25°C

Capacity [Cu. Ft.]	17 Cu. Ft.
Capacity [Liters]	481 Liters
Temperature Range [°C]	-15°C to -25°C
Temperature Range [°F]	5°F to -13°F

Dimensions	Width x Depth x Height
Exterior [Inches]	30.75" x 28.63" x 67.75"
Interior [Inches]	26.13" x 22.5" x 58.25"
Exterior [Centimeters]	78 x 72 x 172
Interior [Centimeters]	66 x 57 x 147

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R290
Exterior	Powder Coating
Interior	Molded ABS
Insulation	Polyurethane Foam
Lockable Door	Keyed
Porthole	3/4" Diameter
Leveling Legs	4
Storage	6 Fixed Shelves
Shipping Weight	235 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 5.5 A
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed Energy Star Rated

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Alarm Relay	Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 2 Year Parts & Labor