

MODEL DHH20-17SDFMS

FLAMMABLE MATERIAL STORAGE FREEZER -15°C TO -25°C

Capacity [Cu. Ft.] 17 Cu. Ft.
Capacity [Liters] 481 Liters
Temperature Range [°C] -15°C to -25°C
Temperature Range [°F] 5°F to -13°F

Dimensions Width x Depth x Height
 Exterior [Inches] 30.75" x 28.5" x 67.75"
 Interior [Inches] 26.25" x 22.5" x 58.25"
 Exterior [Centimeters] 78 x 72 x 172
 Interior [Centimeters] 67 x 57 x 148

Compressor Hermetic Compressor
Air Circulation Gravity Flow
Defrost Method Manual
Refrigerant R290
Exterior Painted Steel
Interior Painted Steel
Insulation High Density Urethane
Lockable Door 1 Door, Keyed Lock
Door Material Insulated Steel
Porthole 3/4" Diameter
Leveling Legs 4
Storage 6 Fixed Shelves
Shipping Weight 220 lbs.

Electrical
Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 6A
NEMA Configuration NEMA 5-15
Certification(s) ETL Listed
 Energy Star Rated

Controller
Display Type Microprocessor, Digital
Alarms High/Low Temp
Min / Max Temperature Display and Reset
Alarm Relay Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» **Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.**

Warranty 2 Years Parts & Labor.