

SO-LOW FLAMMABLE MATERIAL STORAGE FREEZER MODEL DHH20-17SDFMS MANUAL DEFROST

Meets NFPA and OSHA guidelines 45, and 70 for flammable material storage in general laboratory work areas

17 Cubic Ft. Upright -25°C / Solid Door

Specifications:

| | |
|----------------------|---------------------------|
| * Capacity | 17 cubic ft. / 481 liters |
| * Temperature | -15°C to -25°C |
| * Dimensions | |
| * Exterior w x d x h | 30.75" x 28.5" x 67.75" |
| * Interior w x d x h | 26.25" x 22.5" x 58.25" |

Features:

| | |
|----------------------|---|
| *Thermostat | *Microprocessor temperature controller *Remote alarm contacts *Power failure alarm *Audible and visual high/low temperature alarms *Alarm validation (audible, visual, and remote alarm) *Sample simulated digital temperature display |
| *Access Port | 3/8" Access port |
| *Compressor | 1/2 Hp Compressor |
| *Doors | Single solid – key lock |
| *Legs | Four leveling |
| *Shelves | 6 Fixed |
| *Insulation: | High density urethane foam |
| *Refrigerant: | R-600A |
| *Air Circulation | Cold Wall |
| *Defrost | Manual Defrost |
| *Voltage | 115-60-1 |
| *Running Amps | 6 |
| *Supply Plug | 5-15 NEMA |
| *Shipping weight | 220 lbs. |
| *Approvals | ETL Listed |
| *Warranty-(Domestic) | Two year parts and labor. |

