

## MODEL DHH20-20MDP

### MANUAL DEFROST FREEZER -15°C TO -25°C

**Capacity [ Cu. Ft. ]**      **20 Cu. Ft.**  
**Capacity [ Liters ]**      **566 Liters**  
**Temperature Range [ °C ]**      **-15°C to -25°C**  
**Temperature Range [ °F ]**      **5°F to -13°F**

Dimensions      Width x Depth x Height  
 Exterior [ Inches ]      30.5" x 29.5" x 74.75"  
 Interior [ Inches ]      25.5" x 21" x 64.75"  
 Exterior [ Centimeters ]      77 x 74 x 189  
 Interior [ Centimeters ]      64 x 53 x 164

Compressor      Hermetic Compressor  
 Air Circulation      Gravity Flow  
 Defrost Method      Manual  
 Refrigerant      R290  
 Exterior      Powder Coating  
 Interior      Molded ABS  
 Insulation      Polyurethane Foam  
 Lockable Door      Keyed  
 Porthole      3/4" Diameter  
 Leveling Legs      4  
 Storage      6 Fixed Shelves  
 Shipping Weight      275 lbs.

#### Electrical

Voltage:      115 Volts 60 Hertz 1 Phase  
 Amperage Line | Running      15 Amp Dedicated | 5.5 A  
 NEMA Configuration      NEMA 5-15  
 Certification(s)      ETL Listed  
                                          Energy Star Rated

#### Controller

Display Type      Microprocessor, Digital  
 Alarms      High/Low Temp  
 Alarm Relay      Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 2 Year Parts & Labor