

## MODEL DHH20-20SDFMS

# FLAMMABLE MATERIAL STORAGE FREEZER -15°C TO -25°C

**Capacity [ Cu. Ft. ]** 20 Cu. Ft.  
**Capacity [ Liters ]** 566 Liters  
**Temperature Range [ °C ]** -15°C to -25°C  
**Temperature Range [ °F ]** 5°F to -13°F

**Dimensions** Width x Depth x Height  
 Exterior [ Inches ] 30.75" x 29" x 74.5"  
 Interior [ Inches ] 25.75" x 21.75" x 65"  
 Exterior [ Centimeters ] 78 x 74 x 189  
 Interior [ Centimeters ] 65 x 55 x 165

**Compressor** Hermetic Compressor  
**Air Circulation** Gravity Flow  
**Defrost Method** Manual  
**Refrigerant** R290  
**Exterior** Painted Steel  
**Interior** Painted Steel  
**Insulation** High Density Urethane  
**Lockable Door** 1 Door, Keyed Lock  
**Door Material** Insulated Steel  
**Porthole** 3/4" Diameter  
**Leveling Legs** 4  
**Storage** 6 Fixed Shelves  
**Shipping Weight** 220 lbs.

**Electrical**  
**Voltage:** 115 Volts 60 Hertz 1 Phase  
**Amperage Line | Running** 15 Amp Dedicated | 6A  
**NEMA Configuration** NEMA 5-15  
**Certification(s)** ETL Listed  
 Energy Star Rated

**Controller**  
**Display Type** Microprocessor, Digital  
**Alarms** High/Low Temp  
**Min / Max Temperature** Display and Reset  
**Alarm Relay** Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» **Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.**

**Warranty** 2 Years Parts & Labor.

Specifications subject to change without notice