

MODEL DHH20-20SDFX

EXPLOSION PROOF FREEZER -15°C TO -25°C

Capacity [Cu. Ft.]	20 Cu. Ft.
Capacity [Liters]	566 Liters
Temperature Range [°C]	-15°C to -25°C
Temperature Range [°F]	5°F to -13°F

Dimensions	Width x Depth x Height
Exterior [Inches]	30.75" x 29" x 74.5"
Interior [Inches]	25.75" x 21.75" x 65"
Exterior [Centimeters]	78 x 74 x 189
Interior [Centimeters]	65 x 55 x 165

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	High Density Urethane
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Leveling Legs	4
Storage	6 Fixed Shelves
Shipping Weight	240 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6A
Electrical Plug	Not Supplied (Hardwired)
Certification(s)	ETL Listed

Controller

Display Type	Mechanical Thermostat
--------------	-----------------------

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» Meets NFPA / OSHA Guidelines for Hazardous Locations Class 1, Division II, Group C & D

Warranty 1 Year Parts & Labor.