

SO-LOW

Environmental Equipment Co., Inc.

10310 Spartan Dr., Cincinnati, Ohio 45215 USA

(513) 772-9410 e-mail: sales@so-low.com

Fax (513) 772-0570 <http://www.so-low.com>

MODEL DHH20-20SDFX

EXPLOSION PROOF FREEZER -15°C TO -25°C

| | |
|---------------------------------|-----------------------|
| Capacity [Cu. Ft.] | 20 Cu. Ft. |
| Capacity [Liters] | 566 Liters |
| Temperature Range [°C] | -15°C to -25°C |
| Temperature Range [°F] | 5°F to -13°F |

| | |
|--------------------------|------------------------|
| Dimensions | Width x Depth x Height |
| Exterior [Inches] | 30.75" x 29" x 74.5" |
| Interior [Inches] | 25.75" x 21.75" x 65" |
| Exterior [Centimeters] | 78 x 74 x 189 |
| Interior [Centimeters] | 65 x 55 x 165 |

| | |
|-----------------|-----------------------|
| Compressor | Hermetic Compressor |
| Air Circulation | Gravity Flow |
| Defrost Method | Manual |
| Refrigerant | R290 |
| Exterior | Painted Steel |
| Interior | Painted Steel |
| Insulation | High Density Urethane |
| Lockable Door | 1 Door, Keyed Lock |
| Door Material | Insulated Steel |
| Leveling Legs | 4 |
| Storage | 6 Fixed Shelves |
| Shipping Weight | 240 lbs. |

Electrical

| | |
|-------------------------|---------------------------------|
| Voltage: | 115 Volts 60 Hertz 1 Phase |
| Amperage Line Running | 15 Amp Dedicated 6A |
| Electrical Plug | Not Supplied (Hardwired) |
| Certification(s) | ETL Listed Energy Star Rated |

Controller

| | |
|--------------|-----------------------|
| Display Type | Mechanical Thermostat |
|--------------|-----------------------|

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» **Meets NFPA / OSHA Guidelines for Hazardous Locations
Class 1, Division II, Group C & D**

Warranty 1 Year Parts & Labor.

Specifications subject to change without notice