

MODEL DHH20-30SDFMS

FLAMMABLE MATERIAL STORAGE FREEZER -18°C TO -25°C

Capacity [Cu. Ft.]	30 Cu. Ft.
Capacity [Liters]	849 Liters
Temperature Range [°C]	-18°C to -25°C
Temperature Range [°F]	0°F to -13°F

Dimensions	Width x Depth x Height
Exterior [Inches]	35.75 x 35.5 x 72.5"
Interior [Inches]	32" x 22.5" x 65"
Exterior [Centimeters]	90 x 90 x 184
Interior [Centimeters]	81 x 57 x 165

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R404A
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Porthole	3/4" Diameter
Leveling Legs	4
Storage	4 Fixed Shelves
Shipping Weight	380 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed
	Energy Star Rated

Controller

Display Type	Mechanical Thermostat
--------------	-----------------------

Digital Control Option*

*Digital control option includes control & alarm system.
 Note: If a temperature recorder is purchased, the digital control option must also be purchased.

» **Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.**

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & Labor.