

SO-LOW

Environmental Equipment Co., Inc.

10310 Spartan Dr., Cincinnati, Ohio 45215 USA

(513) 772-9410 e-mail: sales@so-low.com

Fax (513) 772-0570 <http://www.so-low.com>

MODEL DHH20-30SDFX

EXPLOSION PROOF FREEZER -18°C TO -25°C

Capacity [Cu. Ft.]	30 Cu. Ft.
Capacity [Liters]	849 Liters
Temperature Range [°C]	-18°C to -25°C
Temperature Range [°F]	0°F to -13°F

Dimensions	Width x Depth x Height
Exterior [Inches]	35.75" x 35" x 77.5"
Interior [Inches]	32" x 22.5" x 65"
Exterior [Centimeters]	90 x 88 x 196
Interior [Centimeters]	81 x 57 x 165

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R404A
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Leveling Legs	4
Storage	4 Fixed Shelves
Shipping Weight	390 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6A
Electrical Plug	Not Supplied (Hardwired)
Certification(s)	UL Listed
	Energy Star Rated

Controller

Display Type	Mechanical Thermostat
Display Housing	Sealed Explosion Proof

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» **Meets NFPA / OSHA Guidelines for Hazardous Locations Class 1, Division II, Group C & D**

Warranty 1 Year Parts & Labor.

Specifications subject to change without notice