

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW FLAMMABLE MATERIAL STORAGE REFRIGERATOR MODEL DHH4-14SDFMS

Meets NFPA and OSHA guidelines 45, and 70 for flammable material storage in general laboratory work areas

14 Cubic Ft. Upright +4°C / Solid Door

Specifications:

* Capacity	14 cubic ft. / 396 liters
* Temperature	1°C to 10°C
* Dimensions	
* Exterior w x d x h	23.5" x 24 ¼" x 73"
* Interior w x d x h	19" x 16.75" x 63.75"

Features:

*Thermostat	*Microprocessor temperature controller *Remote alarm contacts *Power failure alarm *Audible and visual high/low temperature alarms *Alarm validation (audible, visual, and remote alarm) *Sample simulated digital temperature display
*Access Port	3/8" Access port
*Compressor	1/5 Hp Compressor
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	5 adjustable
*Insulation:	High density urethane foam
*Refrigerant:	R-600A
*Air Circulation	Fan
*Defrost	Cycle
*Voltage	115-60-1
*Running Amps	3
*Supply Plug	5-15 NEMA
*Shipping weight	200 lbs.
*Approvals	ETL Listed
*Warranty-(Domestic)	Two year parts and labor.



Specifications subject to change without notice