

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## ***SO-LOW FLAMMABLE MATERIAL STORAGE REFRIGERATOR MODEL DHH4-20SDFMS***

*Meets NFPA and OSHA guidelines 45, and 70 for flammable material storage in general laboratory work areas*

**20 Cubic Ft. Upright +4°C / Solid Door**

### **Specifications:**

* Capacity	20 cubic ft. / 566 liters
* Temperature	1°C to 10°C
* Dimensions	
* Exterior w x d x h	30.75" x 29" x 74.5"
* Interior w x d x h	24.75" x 20.75" x 65.25"

### **Features:**

*Thermostat	*Microprocessor temperature controller
	*Remote alarm contacts
	*Power failure alarm
	*Audible and visual high/low temperature alarms
	*Alarm validation (audible, visual, and remote alarm)
	*Sample simulated digital temperature display
*Access Port	3/8" Access port
*Compressor	1/4 Hp Compressor
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	5 adjustable
*Insulation:	High density urethane foam
*Refrigerant:	R-600A
*Air Circulation	Fan
*Defrost	Cycle
*Voltage	115-60-1
*Running Amps	3
*Supply Plug	5-15 NEMA
*Shipping weight	260 lbs.
*Approvals	ETL Listed
*Warranty-(Domestic)	Two year parts and labor.



Specifications subject to change without notice