

MODEL DHH4-20SDFMS

FLAMMABLE MATERIAL STORAGE REFRIGERATOR 1°C TO 10°C

Capacity [Cu. Ft.] 20 Cu. Ft.
Capacity [Liters] 566 Liters
Temperature Range [°C] 1°C to 10°C
Temperature Range [°F] 33°F to 50°F

Dimensions Width x Depth x Height
 Exterior [Inches] 30.75" x 29" x 74.5"
 Interior [Inches] 24.75" x 20.75" x 65.25"
 Exterior [Centimeters] 78 x 73 x 189
 Interior [Centimeters] 63 x 53 x 166

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Cycle
Refrigerant R600a
Exterior Painted Steel
Interior Painted Steel
Insulation High Density Urethane
Lockable Door 1 Door, Keyed Lock
Door Material Insulated Steel
Porthole 3/4" Diameter
Leveling Legs 4
Storage 5 Adjustable Shelves
Shipping Weight 260 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 3A
NEMA Configuration NEMA 5-15
Certification(s) ETL Listed

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp
Min / Max Temperature Display and Reset
Alarm Relay Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» **Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.**

Warranty 2 Years Parts & Labor.