

MODEL DHK-20RFC

LAB AND PHARMACY REFRIGERATOR FREEZER COMBO

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Total Capacity [Cu. Ft.] 18.33 Cu. Ft. (Total)
Capacity [Liters] 519 Liters (Total)

	REFRIGERATOR	FREEZER
Temperature Range [°C]	2°C to 8°C	-10°C to -26°C
Temperature Range [°F]	35°F to 46°F	14°F to -14.8°F

Dimensions [W x D x H]		
Interior [Inches]	29.5" x 24" x 29"	25.5" x 22" x 22.25"
Interior [Centimeters]	75 x 61 x 74	65 x 56 x 57
Exterior [Inches]	36" x 29.25" x 77.75"	
Exterior [Centimeters]	91 x 74 x 197	


Compressor	Hermetic	Hermetic
Air Circulation	Forced Air	Direct
Defrost Method	Automatic	Manual
Porthole	One 1" Diameter	One 1" Diameter
Refrigerator Storage	3 Adjustable Shelves, Sliding Drawer	
Freezer Storage	3 Fixed Shelves, ABS Drawers	
Refrigerant	Ref: R600a Fzr: R290	
Exterior	Painted Steel	
Interior	Stainless Steel	
Insulation	US EPA and SNAP approved	
Lockable Door	Yes, with separate lock	
Casters	4, with Leveling Feet	
Shipping Weight	500 lbs.	

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 4.55A
Power Consumption	3.76 kW/d
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Type	Microprocessor
Alarms	Two Digital Displays High/Low Temp, Door Ajar Power Failure, Low Battery Sensor Failure, USB Failure
Data Logging	Yes, USB Downloadable, PDF
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts
Backup Battery	Yes, Only Powers Alarm



Warranty 2 Year Parts & Labor