

SO-LOW LABORATORY FREEZER TO -25°C PREMIER MODEL DHK20-20MDP

20.0 Cubic Ft.(566 liters) Upright

Temperature Range: -10°C to -25°C

Solid Door

Specifications:

* Capacity	20.0 cubic ft. / 566 liters
*Temperature Range	-10C to -25°C
* Dimensions	
* Exterior w x d x h	30.5" x 29.5" x 74.75"
* Internal w x d x h	25.5" x 21" x 64.75"

Features:

Temperature Control	Digital display With Hi/Low Alarm
Door	Single solid, locking
Access Port	1/2" Diameter
Interior Finish	Galvanized steel powder coating
Exterior Finish	ABS
Casters	Four leveling legs
Shelves	Six Non-Adjustable
Insulation	Polyurethane foam
Refrigerant:	R600 CFC free
Defrost	Manual Defrost
Voltage: Domestic	115/60/1
Running Amps/Breaker	3.1 / 15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	275 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor

