

## MODEL DHK20-20MDP

### LAB AND PHARMACY FREEZER -15°C TO -25°C

|                                 |                                 |
|---------------------------------|---------------------------------|
| <b>Capacity [ Cu. Ft. ]</b>     | <b>20 Cu. Ft.</b>               |
| <b>Capacity [ Liters ]</b>      | <b>566 Liters</b>               |
| <b>Temperature Range [ °C ]</b> | <b>-15°C to -25°C</b>           |
| <b>Temperature Range [ °F ]</b> | <b>5°F to -13°F</b>             |
| Dimensions                      | Width x Depth x Height          |
| Exterior [ Inches ]             | 30.5" x 29.5" x 74.75"          |
| Interior [ Inches ]             | 25.5" x 21" x 64.75"            |
| Exterior [ Centimeters ]        | 77 x 74 x 189                   |
| Interior [ Centimeters ]        | 64 x 53 x 164                   |
| Compressor                      | Hermetic Compressor             |
| Air Circulation                 | Gravity Flow                    |
| Defrost Method                  | Manual                          |
| Refrigerant                     | R290                            |
| Exterior                        | Steel, Powder Coating           |
| Interior                        | Galvanized Steel, Powder Coated |
| Insulation                      | Polyurethane Foam               |
| Lockable Door                   | Keyed                           |
| Porthole                        | 1/2" Diameter                   |
| Leveling Legs                   | 4                               |
| Storage                         | 6 Fixed Shelves                 |
| Shipping Weight                 | 275 lbs.                        |

#### Electrical

|                         |                            |
|-------------------------|----------------------------|
| Voltage:                | 115 Volts 60 Hertz 1 Phase |
| Amperage Line   Running | 15 Amp Dedicated   3.1A    |
| NEMA Configuration      | NEMA 5-15                  |
| Certification(s)        | ETL Listed                 |

#### Controller

|              |                         |
|--------------|-------------------------|
| Display Type | Microprocessor, Digital |
| Alarms       | High/Low Temp           |
| Alarm Relay  | Yes, Dry Contacts       |

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 2 Year Parts & Labor