

## MODEL DHK4-23GD-T

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

<b>Capacity [ Cu. Ft. ]</b>	<b>23 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>650 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>
<b>Temperature Range [ °F ]</b>	<b>35°F to 46°F</b>
Dimensions	Width x Depth x Height
Exterior [ Inches ]	28.25" x 35" x 78.25"
Interior [ Inches ]	23.75" x 28.5" x 59.75"
Exterior [ Centimeters ]	72 x 89 x 199
Interior [ Centimeters ]	61 x 73 x 152
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Refrigerant	R290
Exterior	Painted Steel
Interior	Stainless Steel
Insulation	US EPA and SNAP approved
Lockable Door	Yes, Keyed
Porthole	Yes, 1" Diameter
Casters	4, with Leveling Feet
Storage	5 Adjustable Shelves,
Shipping Weight	425 lbs.


#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   4.92
Power Consumption	4.54 kW/d
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

#### Controller

Display Type	Microprocessor, Digital
Touchscreen	Yes, Color Display
Alarms	High/Low Temp, Door Ajar Power Failure, Low Battery Sensor Failure, USB Failure
Data Logging	Yes, Adjustable Intervals
Data Download	Yes, via USB, PDF Format
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts
Backup Battery	Yes, Only Powers Alarm

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 2 Year Parts & Labor