

MODEL DHK4-49GD-T

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	46.62 Cu. Ft.
Capacity [Liters]	1320 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	35°F to 46°F
Dimensions	Width x Depth x Height
Exterior [Inches]	57.25" x 34.75" x 78.75"
Interior [Inches]	52.75" x 26.25" x 59.5"
Exterior [Centimeters]	145 x 88 x 200
Interior [Centimeters]	134 x 67 x 152
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Refrigerant	R290
Exterior	Painted Steel
Interior	Stainless Steel
Insulation	US EPA and SNAP approved
Lockable Door	Yes, Keyed
Porthole	Yes, 1" Diameter
Casters	4, with Leveling Feet
Storage	10 Adjustable Shelves
Shipping Weight	800 lbs.


Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 10.06
Power Consumption	11.17 kW/d
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Display Type	Microprocessor, Digital
Touchscreen	Yes, Color Display
Alarms	High/Low Temp, Door Ajar Power Failure, Low Battery Sensor Failure, USB Failure
Data Logging	Yes, Adjustable Intervals
Data Download	Yes, via USB, PDF Format
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts
Backup Battery	Yes, Only Powers Alarm

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 2 Year Parts & Labor