

## MODEL DHK4-49GD

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

**Capacity [ Cu. Ft. ]** 46.62 Cu. Ft.  
**Capacity [ Liters ]** 1320 Liters  
**Temperature Range [ °C ]** 2°C to 8°C  
**Temperature Range [ °F ]** 35°F to 46°F

Dimensions  
Exterior [ Inches ] 57.25" x 34.75" x 78.75"  
Interior [ Inches ] 52.75" x 26.25" x 59.5"  
Exterior [ Centimeters ] 145 x 88 x 200  
Interior [ Centimeters ] 134 x 67 x 152

Compressor Hermetic Compressor  
Air Circulation Forced Air  
Defrost Method Automatic  
Refrigerant R290  
Exterior Painted Steel  
Interior Stainless Steel  
Insulation US EPA and SNAP approved  
Lockable Door Yes, Keyed  
Porthole Yes, 1" Diameter  
Casters 4, with Leveling Feet  
Storage 10 Adjustable Shelves  
Shipping Weight 800 lbs.


#### Electrical

Voltage: 115 Volts 60 Hertz 1 Phase  
Amperage Line | Running 15 Amp Dedicated | 10.06  
Power Consumption 11.17 kW/d  
NEMA Configuration NEMA 5-15  
Certification(s) UL Listed

#### Controller

Display Type Microprocessor, Digital  
Alarms High/Low Temp, Door Ajar  
Power Failure, Low Battery  
Sensor Failure, USB Failure  
Data Logging Yes, Adjustable Intervals  
Data Download Yes, via USB, PDF Format  
Min / Max Temperature Yes, Display and Reset  
Alarm Relay Yes, Dry Contacts  
Backup Battery Yes, Only Powers Alarm

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 2 Year Parts & Labor