

SO-LOW FLAMMABLE MATERIAL STORAGE

FREEZER MODEL DHN20-24FMS

24 Cubic Ft. Upright – Temperature Range: -10°C to -25°C Solid Door

Specifications:

* Capacity	24 cubic ft.
* Temperature Range	-10°C to -25°C
* Dimensions	
* Exterior w x d x h	27.5" x 34.875" x 79.675" 79 x 88 x 202 (cm)
* Interior w x d x h	23.5" x 30" x 59" 59 x 76 x 149 (cm)

Features:

Thermometer	Internal Chamber
Doors	Single Solid, locking
Interior finish	smooth aluminum
Casters	Four
Shelves	Three
Insulation	Polyurethane
Refrigerant:	R404A CFC free
Compressor	3/4 hp
Condenser	Fin and tube
Evaporator	Cold Wall
Thermostat	Adjustable
Defrost	Manual
Voltage: Domestic	115/60/1
Running Amps/Breaker	8.2/15
Supply Plug	5-15P NEMA
Warranty-(Domestic)	UL listed, CSA approved 1 ½ years parts and labor. Compressors additional 3 ½ years, no labor.

