

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW LAB & PHARMACY REFRIGERATOR WITH DRAWERS MODEL DHN4-52GD-D

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

52 Cubic Ft. Upright – Temperature Range: 2°C to 10°C / Glass Doors

Specifications:

* Capacity	52 cubic ft.
* Temperature Range	2°C to 10°C
* Dimensions	
* Exterior w x d x h	55" x 35.5" x 78.625" 139 x 90 x 202 (cm)
* Interior w x d x h	51" x 30" x 59" 129 x 76 x 150 (cm)

Features:

Temperature Control	Digital display with Hi/Low alarm
Doors	Double glass, locking
Interior finish	White
Casters	Four
Drawers	8 stainless steel drawers 18"Wx22"Lx3"H
Interior Lighting	Fluorescent
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/2 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Thermostat	Adjustable
Defrost	Cycle, automatic
Voltage: Domestic	115/60/1
Running Amps/Breaker	10.7/15
Supply Plug	5-15P NEMA UL listed, CSA approved
Shipping Weight	580 lbs.
Warranty-(Domestic)	1 ½ years parts and labor. Compressors additional 3 ½ years, no labor.



Specifications subject to change without notice