

SO-LOW LAB & PHARMACY REFRIGERATOR WITH DRAWERS MODEL DHN4-80GD-D

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

80 Cubic Ft. Upright – Temperature Range: 2°C to 10°C / Glass Doors

Specifications:

*Capacity	80 cubic ft.
*Temperature Range	2°C to 10°C
*Dimensions	
* Exterior w x d x h	82.5” x 35.5” x 79.625” 209 x 90 x 202 (cm)
* Interior w x d x h	78.5” x 30” x 59” 199 x 76 x 150 (cm)

Features:

Temperature Control	Digital display with Hi/Low alarm
Doors	Triple glass, locking
Interior finish	White
Casters	Six
Drawers	12 stainless steel drawers
Interior Lighting	Fluorescent
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	3/4 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Thermostat	Adjustable
Defrost	Cycle, automatic
Voltage: Domestic	115/60/1
Running Amps/Breaker	8.3/15
Supply Plug	L14-20P NEMA
Shipping Weight	UL listed, CSA approved 810 lbs.
Warranty-(Domestic)	1 ½ years parts and labor. Compressors additional 3 ½ years, no labor.

