

MODEL DHP4-23SD

TOUCHSCREEN LAB REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	23 Cu. Ft.
Capacity [Liters]	650 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	36°F to 46°F
Dimensions	Width x Depth x Height
Exterior [Inches]	27" x 32" x 81"
Interior [Inches]	24" x 28.5" x 55"
Exterior [Centimeters]	68 x 81 x 206
Interior [Centimeters]	61 x 72 x 140
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Interior Light	Adjustable LED
Porthole	1/2" Diameter (Top Back Wall)
Casters	4, Locking
Storage	5 Adjustable Shelves
Shipping Weight	440 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 3.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Display Type	Microprocessor, Digital
Touchscreen	Yes, 4.3" Color Display
Alarms	High/Low Temp
Alarm Notification	Yes, Audible, Email, Text
Data Logging	Yes, Adjustable Intervals
Data Download	Yes, via USB, CSV Format
LAN Remote Access	Yes, via VNC
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts

Warranty 3 Years Parts & Labor.
Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

