

MODEL DHP4-27SD

TOUCHSCREEN LAB REFRIGERATOR 2°C TO 8°C

| | |
|---------------------------------|---------------------------------|
| Capacity [Cu. Ft.] | 27 Cu. Ft. |
| Capacity [Liters] | 764 Liters |
| Temperature Range [°C] | 2°C to 8°C |
| Temperature Range [°F] | 36°F to 46°F |
| Dimensions | Width x Depth x Height |
| Exterior [Inches] | 30" x 32" x 81" |
| Exterior [Centimeters] | 76 x 81 x 206 |
| Interior [Inches] | 27" x 28.5" x 55" |
| Interior [Centimeters] | 69 x 72 x 140 |
| Compressor | Hermetic Compressor |
| Air Circulation | Forced Air |
| Defrost Method | Cycle |
| Refrigerant | R290 |
| Exterior | Painted Steel |
| Interior | Painted Steel |
| Insulation | Polyurethane Foam |
| Lockable Door | 1 Door, Keyed Lock |
| Door Material | Insulated Steel |
| Interior Light | Adjustable LED |
| Porthole | 1/2" Diameter (Top Back Wall) |
| Casters | 4, Locking |
| Storage | 5 Adjustable Shelves |
| Shipping Weight | 440 lbs. |

Electrical

| | |
|-------------------------|----------------------------|
| Voltage: | 115 Volts 60 Hertz 1 Phase |
| Amperage Line Running | 15 Amp Dedicated 3.0A |
| NEMA Configuration | NEMA 5-15 |
| Certification(s) | UL Listed |

Controller

| | |
|-----------------------|---------------------------|
| Display Type | Microprocessor, Digital |
| Touchscreen | Yes, 4.3" Color Display |
| Alarms | High/Low Temp |
| Alarm Notification | Yes, Audible, Email, Text |
| Data Logging | Yes, Adjustable Intervals |
| Data Download | Yes, via USB, CSV Format |
| LAN Remote Access | Yes, via VNC |
| Min / Max Temperature | Yes, Display and Reset |
| Alarm Relay | Yes, Dry Contacts |

Warranty 3 Years Parts & Labor.
 Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

