

## MODEL DHP4-38SGD

### TOUCHSCREEN LAB REFRIGERATOR 2°C TO 8°C

<b>Capacity [ Cu. Ft. ]</b>	<b>38 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>1076 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>
<b>Temperature Range [ °F ]</b>	<b>36°F to 46°F</b>
Dimensions	Width x Depth x Height
Exterior [ Inches ]	43.5" x 32" x 81"
Interior [ Inches ]	40.5" x 27" x 55"
Exterior [ Centimeters ]	110 x 81 x 205
Interior [ Centimeters ]	102 x 68 x 140
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	2 Sliding Doors, Keyed Lock
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter
Casters	4, Locking
Storage	10 Adjustable Shelves
Shipping Weight	550 lbs.

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   5.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

#### Controller

Display Type	Microprocessor, Digital
Touchscreen	Yes, 4.3" Color Display
Alarms	High/Low Temp
Alarm Notification	Yes, Audible, Email, Text
Data Logging	Yes, Adjustable Intervals
Data Download	Yes, via USB, CSV Format
LAN Remote Access	Yes, via VNC
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts

**Warranty** 3 Years Parts & Labor.  
Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

