

MODEL DHP4-45SGD

TOUCHSCREEN LAB REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 45 Cu. Ft.
Capacity [Liters] 1286 Liters
Temperature Range [°C] 2°C to 8°C
Temperature Range [°F] 36°F to 46°F

Dimensions Width x Depth x Height
Exterior [Inches] 52" x 32" x 81"
Interior [Inches] 49" x 27" x 55"
Exterior [Centimeters] 132 x 81 x 206
Interior [Centimeters] 125 x 69 x 140

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Cycle
Refrigerant R290
Exterior Painted Steel
Interior Painted Steel
Insulation Polyurethane Foam
Lockable Door 2 Sliding Doors, Keyed Lock
Door Material Glass
Interior Light Adjustable LED
Porthole 1/2" Diameter
Casters 4, Locking
Storage 10 Adjustable Shelves
Shipping Weight 610 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 5.0A
NEMA Configuration NEMA 5-15
Certification(s) UL Listed

Controller

Display Type Microprocessor, Digital
Touchscreen Yes, 4.3" Color Display
Alarms High/Low Temp
Alarm Notification Yes, Audible, Email, Text
Data Logging Yes, Adjustable Intervals
Data Download Yes, via USB, CSV Format
LAN Remote Access Yes, via VNC
Min / Max Temperature Yes, Display and Reset
Alarm Relay Yes, Dry Contacts

Warranty 3 Years Parts & Labor.
Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

