

MODEL DHP4-49GD

TOUCHSCREEN LAB REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 49 Cu. Ft.
Capacity [Liters] 1387 Liters
Temperature Range [°C] 2°C to 8°C
Temperature Range [°F] 36°F to 46°F

Dimensions Width x Depth x Height
Exterior [Inches] 52" x 32" x 81"
Interior [Inches] 47.5" x 28.5" x 55"
Exterior [Centimeters] 132 x 81 x 206
Interior [Centimeters] 120 x 72 x 140

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Cycle
Refrigerant R290
Exterior Painted Steel
Interior Painted Steel
Insulation Polyurethane Foam
Lockable Doors Two Doors, Keyed Lock
Door Material Glass
Interior Light Adjustable LED
Porthole 1/2" Diameter (Top Back Wall)
Casters 4, Locking
Storage 10 Adjustable Shelves
Shipping Weight 610 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 5.0A
NEMA Configuration NEMA 5-15
Certification(s) UL Listed

Controller

Display Type Microprocessor, Digital
Touchscreen Yes, 4.3" Color Display
Alarms High/Low Temp
Alarm Notification Yes, Audible, Email, Text
Data Logging Yes, Adjustable Intervals
Data Download Yes, via USB, CSV Format
LAN Remote Access Yes, via VNC
Min / Max Temperature Yes, Display and Reset
Alarm Relay Yes, Dry Contacts

Warranty 3 Years Parts & Labor.
Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

