

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW LAB & PHARMACY FREEZER

MODEL DHS25-49SD

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

49 Cubic Ft. (1387 liters) Upright

Temperature Range: -18°C to -25°C / 0°F to -13°F
2 Solid Doors

Specifications:

* Capacity	49 cubic ft. / 1387 liters
* Temperature Range	-18°C to -25°C / 0°F to -13° F
* Dimensions	
* Exterior w x d x h	56.69” x 31.5” x 79.25”
* Interior w x d x h	51.97” x 26.18” x 59.84”

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay Battery back up
Access Port	1/2" Diameter
Doors	Two doors solid - lockable
Casters	Heavy duty, swivel
Shelves	Six epoxy coated adjustable
Interior Lights	Led lights
Insulation	Polyurethane CFC Free
Refrigerant:	R-290
Compressor	1 hp
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Automatic
Voltage:	115/60/1
Running Amps/Breaker	6.2/20
Supply Plug	5-20 NEMA
Shipping Weight	450 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice.