

## MODEL DHW20-30

### LAB AND PHARMACY FREEZER 0°C TO -25°C

<b>Capacity [ Cu. Ft. ]</b>	<b>30 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>849 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>0°C to -25°C</b>
<b>Temperature Range [ °F ]</b>	<b>32°F to -13°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	35.5" x 35" x 72"
Interior [ Inches ]	32" x 22.5" x 65"
Exterior [ Centimeters ]	90 x 88 x 182
Interior [ Centimeters ]	81 x 57 x 165

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R404A
Exterior	Painted Steel
Interior	Painted Steel
Insulation	CFC-Free
Lockable Door	Keyed
Door Handle	Stainless Steel
Leveling Legs	4
Storage	4 Fixed Shelves

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   6A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

#### Controller

Display Type	Mechanical Dial
Alarms	None

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 1 Year Parts & Labor