

## MODEL LFB19W1HC

### SOLID DOOR LAB FREEZER -10°C TO -32°C

**Capacity [ Cu. Ft. ] 16.7 Cubic Feet**

**Capacity [ Liters ] 475 Liters**

**Temperature Range [ °C ] -10°C to -32°C**

**Temperature Range [ °F ] +14°F to -26°F**

Dimensions Width x Depth x Height

Exterior [ Inches ] 29.5" x 29.75" x 80.5"

Interior [ Inches ] 24.25" x 19.75" x 55"

Exterior [ Centimeters ] 74 x 75 x 204

Interior [ Centimeters ] 61 x 50 x 139

Compressor Hermetic Compressor

Air Circulation Forced Air

Defrost Method Automatic

Primary Refrigerant(s) R290

Exterior Steel, Painted White

Interior Molded Polystyrol

Insulation Polyurethane Foam

Self-Closing Door 1 Door, Automatic, Recessed

Lockable Door Keyed Lock

Door Material Insulated Steel

Porthole 1" Diameter

Casters 2 Swivel Locking, 2 Fixed

Storage 5 Adjustable Shelves

Shelving Weight Limit 60 Lbs.

Approximate BTU/Hour 1460

Approximate Decibels 62 dB

#### Electrical:

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage Line | Running 15 Amp Dedicated | 5A

NEMA Configuration NEMA 5-15

Certification(s) ETL Listed

#### Controller:

Display Type Microprocessor, Digital

Alarms High / Low Temp

Min / Max Temperature Display and Reset

Alarm Signal Audible and Visual

Alarm Relay Dry Contacts

Alarm Forwarding Yes

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 3 Years Parts & Labor