

MODEL LFBFS05W1HC

UNDERCOUNTER LAB FREEZER -9°C TO -26°C

| | |
|---------------------------------|-----------------------|
| Capacity [Cu. Ft.] | 4.4 Cubic Feet |
| Capacity [Liters] | 123 Liters |
| Temperature Range [°C] | -9°C to -26°C |
| Temperature Range [°F] | +16°F to -15°F |

| | |
|--------------------------|-------------------------|
| Dimensions | Width x Depth x Height |
| Exterior [Inches] | 23.75" x 24.25" x 32.5" |
| Interior [Inches] | 18" x 17.75" x 26.25" |
| Exterior [Centimeters] | 59 x 61 x 82 |
| Interior [Centimeters] | 45 x 45 x 66 |

| | |
|------------------------|-----------------------------|
| Compressor | Hermetic Compressor |
| Air Circulation | Gravity Flow |
| Defrost Method | Manual |
| Primary Refrigerant(s) | R600 |
| Exterior | Steel, Painted White |
| Interior | Molded Polystyrol |
| Insulation | Polyurethane Foam |
| Self-Closing Door | 1 Door, Automatic, Recessed |
| Lockable Door | Keyed Lock |
| Door Material | Insulated Steel |
| Leveling Legs | 4 |
| Shelving | 3 Fixed Drawers, 1 Basket |
| Shelving Weight Limit | 50 Lbs. |
| Approximate BTU/Hour | 730 |
| Approximate Decibels | 45 dB |

Electrical:

| | |
|-------------------------|----------------------------|
| Voltage: | 115 Volts 60 Hertz 1 Phase |
| Amperage Line Running | 15 Amp Dedicated 2.5A |
| NEMA Configuration | NEMA 5-15 |
| Certification(s) | ETL Listed |

Controller:

| | |
|-----------------------|-------------------------|
| Temperature Control | Microprocessor, Digital |
| Alarms | High / Low Temp |
| Min / Max Temperature | Display and Reset |
| Alarm Signal | Audible and Visual |
| Alarm Relay | Dry Contacts |
| Alarm Forwarding | Yes |

Interior is free from sources of ignition and meets the requirements for Class I, Division 2, for storage of flammable substances in closed containers. Product complies with the application requirements of NFPA 45 and NFPA 99.

Warranty 3 Years Parts & Labor

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

