

## MODEL LFT30W1HC

### SOLID DOOR LAB FREEZER -10°C TO -35°C

<b>Capacity [ Cu. Ft. ]</b>	<b>29.1 Cubic Feet</b>
<b>Capacity [ Liters ]</b>	<b>831 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>-10°C to -35°C</b>
<b>Temperature Range [ °F ]</b>	<b>+14°F to -31°F</b>
Dimensions	Width x Depth x Height
Exterior [ Inches ]	31.25" x 38.75" x 85.25"
Interior [ Inches ]	24.75" x 31.75" x 59.25"
Exterior [ Centimeters ]	79 x 98 x 216
Interior [ Centimeters ]	62 x 80 x 150
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Primary Refrigerant(s)	R290
Exterior	Steel, Painted White
Interior	Stainless Steel
Insulation	Polyurethane Foam
Interior Lighting	LED Ceiling Light
Self-Closing Door	1 Door, Automatic, Recessed
Lockable Door	Keyed Lock
Door Material	Insulated Steel
Porthole	1" Diameter
Casters	2 Swivel Locking, 2 Fixed
Shelving	4 Adjustable Shelves
Shelving Weight Limit	132 Lbs.
Approximate BTU/Hour	2930
Approximate Decibels	65 dB

#### Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   10A
NEMA Configuration	NEMA 5-15
Certification(s)	CSA Listed Energy Star Rated

#### Control:

Display Type	Microprocessor, Digital
Alarms	High / Low Temp
Min / Max Temperature	Display and Reset
Alarm Signal	Audible and Visual
Alarm Relay	Dry Contacts
Alarm Forwarding	Yes

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 3 Years Parts & Labor