

MODEL LFT50W2HC

DOUBLE SOLID DOOR LAB FREEZER -10°C TO -26°C

Capacity [Cu. Ft.]	46.4 Cubic Feet
Capacity [Liters]	1325 Liters
Temperature Range [°C]	-10°C to -26°C
Temperature Range [°F]	+14°F to -15°F

Dimensions	Width x Depth x Height
Exterior [Inches]	56.5" x 32.75" x 85"
Interior [Inches]	48.5" x 25.5" x 59"
Exterior [Centimeters]	143 x 83 x 216
Interior [Centimeters]	123 x 65 x 150

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Primary Refrigerant(s)	R290
Exterior	Steel, Painted White
Interior	Stainless Steel
Insulation	Polyurethane Foam
Interior Lighting	LED Ceiling Light
Self-Closing Door	2 Doors, Automatic, Recessed
Lockable Door	Keyed Lock
Door Material	Insulated Steel
Porthole	1" Diameter
Casters	3 Swivel Locking, 3 Fixed
Storage	8 Adjustable Shelves
Shelving Weight Limit	132 Lbs.
Approximate BTU/Hour	2930
Approximate Decibels	65 dB

Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 10A
NEMA Configuration	NEMA 5-15
Certification(s)	CSA Listed
	Energy Star Listed

Control:

Temperature Control	Microprocessor, Digital
Alarms	High / Low Temp
Min / Max Temperature	Display and Reset
Alarm Signal	Audible and Visual
Alarm Relay	Dry Contacts
Alarm Forwarding	Yes

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 3 Years Parts & Labor