

MODEL LRB05G1HC

UNDERCOUNTER LAB REFRIGERATOR 3°C TO 16°C

Capacity [Cu. Ft.] 5.1 Cubic Feet
Capacity [Liters] 143 Liters
Temperature Range [°C] 3°C to 16°C
Temperature Range [°F] 37°F to 61°F

Dimensions Width x Depth x Height
Exterior [Inches] 23.75" x 24.5" x 32.5"
Interior [Inches] 17.25" x 17.25" x 26.25
Exterior [Centimeters] 59 x 61 x 82
Interior [Centimeters] 44 x 43 x 67

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Automatic
Primary Refrigerant(s) R600a
Exterior Steel, Painted White
Interior Molded Polystyrol
Insulation Polyurethane Foam
Interior Light LED Ceiling Light
Self-Closing Door 1 Door, Automatic, Recessed
Lockable Door Keyed Lock
Door Material Glass
Porthole 1/2" Diameter
Leveling Legs 4
Storage 3 Adjustable Shelves
Shelving Weight Limit 45 Lbs.
Approximate BTU/Hour 440
Approximate Decibels 52 dB

Electrical:

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 1.5A
NEMA Configuration NEMA 5-15
Certification(s) ETL Listed

Control:

Temperature Control Microprocessor, Digital
Alarms High / Low Temp
Min / Max Temperature Display and Reset
Alarm Signal Audible and Visual
Alarm Relay Dry Contacts
Alarm Forwarding Yes

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 3 Years Parts & Labor