

MODEL LRB05W1HC

UNDERCOUNTER LAB REFRIGERATOR 3°C TO 16°C

Capacity [Cu. Ft.]	4.7 Cubic Feet
Capacity [Liters]	133 Liters
Temperature Range [°C]	3°C to 16°C
Temperature Range [°F]	37°F to 61°F

Dimensions	Width x Depth x Height
Exterior [Inches]	23.75" x 24.5" x 32.5"
Interior [Inches]	17.25" x 17.25" x 26.25
Exterior [Centimeters]	59 x 61 x 82
Interior [Centimeters]	44 x 43 x 67
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Primary Refrigerant(s)	R600a
Exterior	Steel, Painted White
Interior	Molded Polystyrol
Insulation	Polyurethane Foam
Interior Light	LED Ceiling Light
Self-Closing Door	1 Door, Automatic, Recessed
Lockable Door	Keyed Lock
Door Material	Insulated Steel
Porthole	1/2" Diameter
Leveling Legs	4
Storage	3 Adjustable Shelves
Shelving Weight Limit	45 Lbs.
Approximate BTU/Hour	440
Approximate Decibels	52 dB

Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 1.5A
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed

Control:

Temperature Control	Microprocessor, Digital
Alarms	High / Low Temp
Min / Max Temperature	Display and Reset
Alarm Signal	Audible and Visual
Alarm Relay	Dry Contacts
Alarm Forwarding	Yes

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 3 Years Parts & Labor