

MODEL LRB23W1HC

SOLID DOOR LAB REFRIGERATOR 3°C TO 16°C

Capacity [Cu. Ft.]	22 Cubic Feet
Capacity [Liters]	627 Liters
Temperature Range [°C]	3°C to 16°C
Temperature Range [°F]	37°F to 61°F
Dimensions	Width x Depth x Height
Exterior [Inches]	29.5" x 29.75" x 80.5"
Interior [Inches]	25.5" x 22.75" x 67.25"
Exterior [Centimeters]	74 x 75 x 204
Interior [Centimeters]	65 x 58 x 171
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Primary Refrigerant(s)	R290
Exterior	Steel, Painted White
Interior	Molded Polystyrol
Insulation	Polyurethane Foam
Interior Light	LED Light
Self-Closing Door	1 Door, Automatic, Recessed
Lockable Door	Keyed Lock
Door Material	Insulated Steel
Porthole	1/2" Diameter
Support	2 Swivel Locking, 2 Fixed
Storage	6 Adjustable Shelves
Shelving Weight Limit	60 Lbs.
Approximate BTU/Hour	880
Approximate Decibels	55 dB

Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 3A
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed

Control:

Temperature Control	Microprocessor, Digital
Alarms	High / Low Temp
Min / Max Temperature	Display and Reset
Alarm Signal	Audible and Visual
Alarm Relay	Dry Contacts
Alarm Forwarding	Yes

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 3 Years Parts & Labor