


LRBFS06W1HC



Unit Dimensions



Interior is free from sources of ignition and meets the requirements for Class I, Division 2 for storage of flammable substances in closed containers. Product complies with the applicable requirements of NFPA 45 and NFPA 99.

Subject to alterations.

Laboratory refrigerator with spark-free interior	LRBFS06W1HC
Total gross/net capacity l (cu.ft.)	180/165 (6.4/5.8)
Exterior dimensions in mm (w/d/h)	597/611/859
in inches (w/d/h)	23 9/16"/24 1/16"/33 7/8"
Interior dimensions in mm (w/d/h)	513/441/702
in inches (w/d/h)	20 1/4"/17 3/8"/27 11/16"
Energy consumption in 365 days *	372 kWh
Ambient temperature in °C (°F)	10°C to 27°C (50°F to 81°F)
Refrigerant	R600a
Sound rating **	50 dB(A)
Connection rating	1.5 A
Frequency/voltage	60 Hz/115V~
Cooling system	forced-air
Defrost	automatic
Temperature range in °C (°F)	1°C to 15°C (34°F to 59°F)
External cabinet finish	steel/white
Door/lid material	steel
Internal liner material	commercial grade molded polystyrol
Type of control	analog control
Temperature display	external digital
Interior light	n/a
Adjustable shelves	3
Shelf material	glass
Shelf load weight in kg (lbs)	40 (88.2)
Usable shelf area in mm (inches)	513/412 (20 1/4"/16 1/4")
Handle	ergonomic slimline handle
Lock	recessed
Door hinges	right, reversible on site
Insulation in mm (inches)	16–55 (11/16"–2 3/16")
Gross/net weight in kg (lbs)	42/39 (92.6/86)
Accessories	
Glass shelf	9293599

* Measured unloaded at 26.6°C (80°F) ambient temperature with set point of 3.3°C (38°F).

** Noise level – sound power re 1pW – according to IEC 60704-2-14.

SO-LOW

AUTHORIZED U.S. DEALER OF
LIEBHERR

WWW.SO-LOW.COM • SALES@SO-LOW.COM • 10310 SPARTAN DRIVE, CINCINNATI, OH 45215 • 513-772-9410