

MODEL LRT50G2HC

TWO GLASS DOOR LAB REFRIGERATOR 2°C TO 16°C

Capacity [Cu. Ft.]	47.5 Cubic Feet
Capacity [Liters]	1357 Liters
Temperature Range [°C]	2°C to 16°C
Temperature Range [°F]	35°F to 61°F

Dimensions	Width x Depth x Height
Exterior [Inches]	56.5" x 32.75" x 85.25"
Interior [Inches]	48.5" x 25.5" x 59"
Exterior [Centimeters]	143 x 83 x 216
Interior [Centimeters]	126 x 65 x 150

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Primary Refrigerant(s)	R290
Exterior	Steel, Painted White
Interior	Stainless Steel
Insulation	Polyurethane Foam
Interior Lighting	LED Ceiling Light
Self-Closing Door	2 Doors, Automatic, Recessed
Lockable Door	Keyed Lock
Door Material	Glass
Porthole	1/2" Diameter
Casters	3 Swivel Locking, 3 Fixed
Storage	8 Adjustable Shelves
Shelving Weight Limit	60 Lbs.
Approximate BTU/Hour	1465
Approximate Decibels	62 dB

Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 5A
NEMA Configuration	NEMA 5-15
Certification(s)	CSA Listed Energy Star Rated

Control:

Temperature Control	Microprocessor, Digital
Alarms	High / Low Temp
Min / Max Temperature	Display and Reset
Alarm Signal	Audible and Visual
Alarm Relay	Dry Contacts
Alarm Forwarding	Yes

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 3 Years Parts & Labor