

## MODEL LRT50W2HC

### TWO SOLID DOOR LAB REFRIGERATOR 2°C TO 16°C

<b>Capacity [ Cu. Ft. ]</b>	<b>46.4 Cubic Feet</b>
<b>Capacity [ Liters ]</b>	<b>1325 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 16°C</b>
<b>Temperature Range [ °F ]</b>	<b>35°F to 61°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	56.75" x 32.7" x 85.25"
Interior [ Inches ]	48.5" x 25.5" x 59"
Exterior [ Centimeters ]	143 x 83 x 216
Interior [ Centimeters ]	126 x 65 x 150

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Primary Refrigerant(s)	R290
Exterior	Steel, Painted White
Interior	Stainless Steel
Insulation	Polyurethane Foam
Interior Lighting	LED Ceiling Light
Self-Closing Door	2 Doors, Automatic, Recessed
Lockable Door	Keyed Lock
Door Material	Insulated Steel
Porthole	1/2" Diameter
Casters	2 Swivel Locking, 2 Fixed
Shelving	8 Adjustable Shelves
Shelving Weight Limit	60 Lbs.
Approximate BTU/Hour	1465
Approximate Decibels	62 dB

#### Electrical:

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   5A
NEMA Configuration	NEMA 5-15
Certification(s)	CSA Listed Energy Star Rated

#### Control:

Temperature Control	Microprocessor, Digital
Alarms	High / Low Temp
Min / Max Temperature	Display and Reset
Alarm Signal	Audible and Visual
Alarm Relay	Dry Contacts
Alarm Forwarding	Yes

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 3 Years Parts & Labor