

## MODEL MV20-1.5UCF

### UNDERCOUNTER FREEZER TO -20°C

<b>Capacity [ Cu. Ft. ]</b>	<b>1.5 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>42 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>-20°C</b>
<b>Temperature Range [ °F ]</b>	<b>-4°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	18.5" x 19.5" x 19.5"
Interior [ Inches ]	14" x 13.25" x 14.5"
Exterior [ Centimeters ]	46 x 49 x 49
Interior [ Centimeters ]	35 x 33 x 36

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R600a
Exterior	Painted Steel
Interior	Steel
Insulation	Urethane Foam
Lockable Door	Keyed
Porthole	Yes, 3/8" Diameter
Door Handle	Stainless Steel
Leveling Legs	4
Storage	1 Fixed Shelf

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   1.1A
NEMA Configuration	NEMA 5-15

#### Controller

Display Type	Mechanical Dial
--------------	-----------------



**Warranty** 1 Year Parts & Labor