

SO-LOW UNDERCOUNTER FREEZER MODEL MV20-1.5UCF

**1.5 Cubic Ft. Upright – Temperature Range: -20°C
One Solid Door**

Specifications:

* Capacity	1.5 cubic ft.
* Temperature Range	-20°C
* Dimensions	
* Exterior w x d x h	18.5" x 19.5" x 19.5" 46 x 49 x 49 (cm)
* Interior w x d x h	14" x 13.25" x 14.5" 35 x 33 x 36 (cm)

Features:

Standard undercounter freezer (freestanding)
Mechanical thermostat
Average cycle temperature: -20°C
Keyed door lock
Magnetic door gasket for positive seal
Leveling legs
Probe access port (3/8")
White cabinet and door
High density urethane foam
CFC-Free refrigerant (R134a)
Manual defrost
115V, 60 Hz, 1/10 HP, 1.1 Amps, 310 BTUs/hr
Shipping weight: 50 lbs
Two year parts and labor warranty,
plus an additional three year parts warranty on compressor.

