

MODEL MV20-2UCF

LAB AND PHARMACY FREEZER -10°C TO -25°C

| | |
|---------------------------------|----------------------------------|
| Capacity [Cu. Ft.] | 1.7 Cu. Ft. |
| Capacity [Liters] | 48 Liters |
| Temperature Range [°C] | -15°C to -25°C |
| Temperature Range [°F] | 5°F to -13°F |
| Dimensions | Width x Depth x Height |
| Exterior [Inches] | 28.75" x 24" x 20.5" |
| Interior [Inches] | 17.75" x 13.75" x 16" |
| Exterior [Centimeters] | 74 x 61 x 53 |
| Interior [Centimeters] | 45 x 34 x 40 |
| Compressor | Hermetic Compressor |
| Air Circulation | Direct Cooling |
| Defrost Method | Manual |
| Refrigerant | R600a |
| Exterior | Painted Steel |
| Interior | Painted Aluminum |
| Insulation | US EPA and SNAP approved |
| Lockable Door | Yes, Keyed |
| Porthole | Yes, 3/8" Diameter |
| Casters | None |
| Leveling Legs | Yes, 4 |
| Storage | 2 Adjustable Shelves, Base Shelf |
| Shipping Weight | 115 lbs. |


Electrical

| | |
|-------------------------|------------------------------|
| Voltage: | 115 Volts 60 Hertz 1 Phase |
| Amperage Line Running | 15 Amp Dedicated 1A |
| NEMA Configuration | NEMA 5-15 |
| Certification(s) | UL Listed, C-UL, Energy Star |

Controller

| | |
|-----------------------|---|
| Display Type | Microprocessor, Digital |
| Alarms | High/Low Temp, Door Ajar Power Failure, Low Battery Sensor Failure, USB Failure |
| Data Logging | Yes, Adjustable Intervals |
| Data Download | N/A |
| Min / Max Temperature | Yes, Display and Reset |
| Alarm Relay | Yes, Dry Contacts |
| Backup Battery | N/A |

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 2 Year Parts & Labor