

SO-LOW UNDERCOUNTER FREEZER FLAMMABLE MATERIAL STORAGE MODEL MV23-4UCFMSF

4.0 Cubic Ft.(114 liters) Upright

Temperature Range: -15°C to -25°C

Solid Door

Specifications:

- * Capacity 4.0 cubic ft. / 114 liters
- * Temperature Range -15°C to -25°C
- * Dimensions
- * Exterior w x d x h 23-5/8" x 24-1/2" x 32-3/16"

Features:

- *DIGITAL TEMPERATURE DISPLAY
- *AUDIBLE AND VISUAL HIGH AND LOW ALARM
- *ADDITIONAL 1" CLEARANCE FROM TOP AND 4" FROM BACK FOR PROPER VENTILATION REQUIRED WHEN RECESSED
- * HFC-FREE R600 REFRIGERATION SYSTEM
- *ACCESS PORT 3/8"
- *115 VOLT, 60 HERTZ, 1 PHASE
- *1.4 AMPS
- *LISTED IN UNDERWRITERS LABORATORIES FILE SA7488(N)
- *MEETS NATIONAL FIRE PROTECTION ASSOCIATION STANDARD 99
- *COMPLIES WITH OSHA ARTICLE NO. CFR1910.308
- *ONE YEAR PARTS AND LABOR WARRANTY ADDITIONAL FOUR YEARS ON COMPRESSORS



Specifications subject to change without notice