

## *SO-LOW UNDERCOUNTER FREEZER MODEL MV30-4UCF*

**4.0 Cubic Ft. Upright – Temperature Range: 0°C to -30°C**

**One Solid Door**

***Specifications:***

* Capacity	4.0 cubic ft.
* Temperature Range	0°C to -30°C
* Dimensions	
* Exterior w x d x h	23.75" x 21.5" x 32" 60 x 64 x 84 (cm)
* Interior w x d x h	19.5" x 15.75" x 22" 50 x 40 x 56 (cm)

***Features:***

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts
Access Port	1/2"
Door	Single solid, locking
Shelves	Two adjustable, 1 full, 1 half
Insulation	Polyurethane foam
Refrigerant:	R600
Compressor	1/3 hp
Defrost	Manual
Voltage: Domestic	115/60/1
Running Amps/Breaker	2/15
Supply Plug	5-15P NEMA
Shipping Weight	150 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice