

MODEL MV4-6UCFMSR

FLAMMABLE MATERIAL STORAGE REFRIGERATOR 1°C TO 10°C

Capacity [Cu. Ft.]	5 Cu. Ft.
Capacity [Liters]	414 Liters
Temperature Range [°C]	1°C to 10°C
Temperature Range [°F]	33°F to 50°F

Dimensions	Width x Depth x Height
Exterior [Inches]	23.75" x 24.25" x 32.25"
Interior [Inches]	20.25" x 17" x 27"
Exterior [Centimeters]	60 x 61 x 81
Interior [Centimeters]	51 x 43 x 68

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Cycle
Refrigerant	R600a
Exterior	Painted Steel
Interior	Painted Steel
Insulation	High Density Urethane
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Porthole	3/8" Diameter
Leveling Legs	4
Storage	3 Adjustable Shelves
Shipping Weight	120 lbs.
Clearance Requirements	4" Clearance on Sides & Back

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 1.2A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed (SA7488(N))

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts
Backup Battery	Yes, Only Powers Alarm

- » **Meets National Fire Protection Assoc. Standard 99**
- » **Complies with OSHA Article No. CFR1910.308**
- » **Enclosed Overload, Relay, and Thermostat**
- » **Static Free Electric Condenser Fan**

Warranty 2 Years Parts & Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

