

SO-LOW UNDERCOUNTER REFRIGERATOR EXPLOSION PROOF - MODEL MV4-6UCXR

6.1 Cubic Ft.(173 liters) Upright

Temperature Range: 1°C to 7°C

Solid Door

Specifications:

* Capacity	6.1 cubic ft. / 173 liters
* Temperature Range	1°C to 7°C
* Dimensions	
* Exterior w x d x h	23-7/8" x 26" x 34" 61 x 66 x 86 (cm)

Features:

- *MANUAL DEFROST
- *STATIC CONDENSER
- *ADDITIONAL 1" CLEARANCE FROM TOP
AND 1.5" FROM BACK REQUIRED FOR
PROPER VENTILATION WHEN RECESSED
- *R134A REFRIGERATION SYSTEM
- *3 IN-DOOR SHELVES
- *115 VOLT, 60 HERTZ, 1 PHASE
- *2.5 AMPS
- *C-UL LISTED

Specifications subject to change without notice



Explosion proof units are approved for use in Class 1, Group C and D environments, NFPA Standard 497A