

## *SO-LOW UNDERCOUNTER FREEZER MODEL MV85-2*

**2 Cubic Ft. (54 liters) Upright – Temperature Range: -40°C to -85°C  
Solid Door**

### ***Specifications:***

* Capacity	2 cubic ft. / 54 liters
* Temperature Range	-40°C to -85°C
* Dimensions	
* Exterior w x d x h	23" x 27" x 32" 60 x 70 x 81 (cm)
* Interior w x d x h	14" x 18" x 12.5" 36 x 46 x 32 (cm)


### ***Features:***

Temperature Control	Digital – Battery Backup Datalogger – Hi/Low Alarm
Clearance	4" on sides and back
Doors	Single Solid, locking
Legs	Four, leveling
Insulation	Polyurethane CFC Free
Evaporator	Cold wall
Defrost	Manual
Voltage: Domestic	115/60/1 – 6.3 amps
Supply Plug	NEMA 5-15P ETL listed
Shipping Weight	265 lbs.
Warranty	Two year parts and labor

### ***Inventory Racks:***

#### **Rack # 17-15-2**

Rack holds fifteen  
2" high x 5.25" x 5.25"  
boxes with dividers.  
MV85-2 holds 2 racks.



Specifications subject to change without notice