

MVK4-2UCRGDDA

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	2.5 Cu. Ft.
Capacity [Liters]	70 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	35°F to 46°F

Dimensions	Width x Depth x Height
Exterior [Inches]	20" x 21.5" x 25"
Interior [Inches]	17" x 15" x 17"
Exterior [Centimeters]	50 x 54 x 63
Interior [Centimeters]	43 x 38 x 43

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Automatic
Refrigerant	R600a
Exterior	Painted Steel, Glass Door
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	Keyed
Porthole	1/2" Diameter
Leveling Feet	4
Storage	2 Adjustable Shelves
Shipping Weight	110 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 1.2A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temperature
Alarm Relay	Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & Labor